VR8V
VALVE

®

CORA

®

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







The VR8V Valve is an aseptic valve, designed and manufactured to dose powder and granules in the presence of pressure and vacuum. The VR8V Valve allows continuous powder and granule flow thus obtaining an extremely uniform dosing during the feeding of Mills, Press Machines, Filling Machines, Micro dosing Machines and Sieving Machines. The CO.RA. VR8V valve has a spherical design that allows to intercept the charging system with the discharging line of the product.

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For further information please refer to our "Pharmaceutical Process" brochure

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For further information please refer to our "Chemical Process" brochure

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TECHNICAL FEATURES


- Half-bodies in AISI 316L (DIN 1.4404) satin finish (mirror polished upon request)
- Rotor in AISI 316L (DIN 1.4404) mirror polished
- 8 pocket rotor with variable speed to dose a controlled amount of product
- NO O-ring on the shaft
- Easy to disassemble, inspect and clean very light weight and compact
- Materials: stainless steel and elastomers (FDA approved, certificates on request)
- Manufactured according to G.M.P.
- It allows continuous feeding of granules and powders without obstructing the flow
- Traceability in time (batch number)
- Compact and light weight compared to conventional rotary valves
- The valve always prevents dust leakage non matter the position of the rotor

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Dimensions and technical details can be subject to variation without notice.

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H4, H5 and D2 may vary depending on choice of motor.

Diametro valvola Valve diameter	Volume lobo dm ³ Pocket volume dm ³	Volume di un giro completo dm ³ Volume of a full revolution dm ³	Capacità massima di dosaggio dm ³ 60 giri/minuto (massima velocità) Maximum dosing capacity dm ³ 60 revolutions / minute (maximum speed)
80	0,023	0,184	662,4
100	0,053	0,424	1.526,4
* 150	0,171	1,368	4.924,8

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The motors for the VR8V and VR8VG can be: pneumatic, electric brushless motor with stainless steel body or stainless steel casing and stepper motor.

If the valve is motorized with a pneumatic motor, it is not possible to control the number of revolutions, so there is no possibility of obtaining a controlled flow rate. In such case you would only have a continuous flow of product, the throughput depends on the size of the valve and on the number of revolutions calculated in the design phase.

By using a brushless motor, stepper motor or a three-phase electric motor powered by a variable frequency drive (inverter) it is possible to perform high precision product dosing, since it is possible to check, and then vary the number of revolutions of the valve. The flow of product through the valve can vary from maximum to minimal flow (see table above).

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The dosing accuracy of the VR8V and VR8VG can get to ± 5 grams with 90% of the pharmaceutical products, to achieve this result CO.RA. continuously improves and develops its products and has elaborated the valve rotor in order to improve the dosing precision of the DN 80 and DN 100, while for the DN 150 the flow rate of the rotor has been increased in order to provide the dosing of significant quantities of product in a short period of time. These high-precision results can be achieved with the above mentioned rotors, but only if used with a motor and control and weighing system supplied by CO.RA. CO.RA. has studied a solution to achieve an even higher dosing accuracy and it is to customize the volume and the number of pockets of the valve rotor which allow it to reach a dosing accuracy of $\pm 2 - 3$ grams.

In these cases, the valve requires fine tuning based on dosing tests carried out with the customer's product. These tests can be performed with the VR8V and VR8VG available in our Technology Center. The possibility of customizing the rotor both in volume in number of pockets ensures great performance flexibility, allowing the optimization of the process at any time simply by replacing the rotor.

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THE VR8VG AND ITS GASKETS

The VR8VG is a variant of VR8V and it is designed to guarantee a superior seal thanks to the sealing gasket between the half-bodies. The gasket types currently available for this valve are: P.T.F.E. core with white silicone, P.T.F.E. core with black silicone.

The VR8VG is designed to ensure a better seal if it is controlled by a brushless motor which positions the pockets in axis with the gasket during the arresting phase. The half-body version and drives are the same as for the VR8V, as are the overall dimensions and throughputs. La VR8VG is produced upon request. The seals are molded in a cleanroom and respect the current FDA and GMP regulations.

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DOCUMENTATION AND CERTIFICATES

The documentation supplied, upon request, with the valve is:

- FDA, USP and 3.1 material certificates of the parts in contact with product
- Use and Maintenance Manual
- Certificate of surface roughness
- Certificate of conformity CE
- FAT and SAT protocols
- IQ and OQ Protocol
- Other upon request

The VR8V and VR8VG complies with Machinery Directive, is CE marked and also can be supplied for Atex Zone. The Atex marking area is following:

II 1GD/2GD c 85°C (T6) 0°C Ta+40°
 II 1D/2GD c 85°C (T6) 0°C Ta+40°
 II 2GD/2GD c IIA 85°C (T6) 0°C Ta+40°
 II 2GD/2GD c IIB 85°C (T6) 0°C Ta+40°

Services related to the supply

The services that can be related to the supply are:

- Assistance within 48 hours from malfunctioning

The VR8V and VR8VG is designed and manufactured respecting customer requirements, and to better guarantee this service it is possible to make tests at our premises with one of our VR8V and VR8VG in our Technology Center.

USEFUL INFORMATION

For further information please contact our Sales Dept by phone +39 0583 20590 or by e-mail: sales@coraitaly.net

UNI EN ISO 9001:2008 - Cert. No. 5762/1

REG. No. IT-71044

ATEX 94/09/CE

CO.RA. s.r.l. - 55011 Spianate - Altopascio (Lucca) - Italy - Loc. Chiappini, 51
Tel. +39 0583 20590 r.a. - Fax +39 0583 20481 - info@coraitaly.net - www.coraitaly.net